

Anti-Recombinant-HA-Tag antibody [SR0595] (STJA0010595)

STJA0010595

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit monoclonal antibody anti-Recombinant-HA-Tag is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	HA-Tag

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SR0595
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB (Species independent) 1:1000-1:2000
Formulation	PBS with 0.15% ProClin300, 100 Mu g/mL BSA and 50% glycerol.
Isotype	IgG
Storage Instruction	Store at -20C for up to one year, and avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	KLH conjugated Synthetic peptide corresponding to YPYDVPDYA
Immunogen Region	
Specificity	
Immunogen Sequence	

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081