

Anti-Recombinant-mCherry antibody [SM0626] (STJA0010626)

STJA0010626

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Recombinant-mCherry is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Mouse
Reactivity	mCherry-Tag

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SM0626
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB (Species independent) 1:500-1:2000
Formulation	PBS with 0.15% ProClin300, 100 µg/mL BSA and 50% glycerol.
Isotype	IgG1
Storage Instruction	Store at -20°C for up to one year, and avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Recombinant protein corresponding to mCherry
Immunogen Region	
Specificity	
Immunogen Sequence	