

Anti-Recombinant-CDK2 antibody [SM0628] (STJA0010628)

STJA0010628

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Recombinant-CDK2 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SM0628
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB (H, M, R) 1:500-1:1000
Formulation	PBS with 0.15% ProClin300, 100 µg/mL BSA and 50% glycerol.
Isotype	IgG1
Storage Instruction	Store at -20°C for up to one year, and avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Gene ID	12566
Gene Symbol	Cdk2
Uniprot ID	CDK2_MOUSE
Immunogen	KLH conjugated Synthetic peptide corresponding to Mouse Cdk2
Immunogen Region	
Specificity	
Immunogen Sequence	