

Anti-Recombinant-Myeloperoxidase antibody [SM0631] (STJA0010631)

STJA0010631

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Recombinant-Myeloperoxidase is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
Applications	IHC/IF
Host/Source	Mouse
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SM0631
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC/IF (M, R) 1:500-1:1000
Formulation	PBS with 0.15% ProClin300, 100 Mu g/mL BSA and 50% glycerol.
Isotype	IgG1
Storage Instruction	Store at -20C for up to one year, and avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Gene ID	17523
Gene Symbol	Mpo
Uniprot ID	PERM_MOUSE
Immunogen	KLH conjugated Synthetic peptide corresponding to Mouse MPO
Immunogen Region	
Specificity	
Immunogen Sequence	