

Anti-Recombinant-COX IV antibody [SM0645] (STJA0010645)

STJA0010645

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Recombinant-COX IV is suitable for use in Western Blot, Immunohistochemistry and Immunofluorescence research applications.
Applications	WB/IHC/IF
Host/Source	Mouse
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SM0645
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB (M, R) 1:1000-1:2000 IHC/IF (M, R) 1:500-1:1000
Formulation	PBS with 0.15% ProClin300, 100 µg/mL BSA and 50% glycerol.
Isotype	IgG1
Storage	Store at -20°C for up to one year, and avoid repeated freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	12857
Gene Symbol	Cox4i1
Uniprot ID	COX41_MOUSE
Immunogen	Recombinant protein corresponding to Mouse COX IV
Immunogen Region	
Specificity	
Immunogen Sequence	