

Anti-Recombinant-Claudin 1 antibody [SR0646] (STJA0010646)

STJA0010646

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit monoclonal antibody anti-Recombinant-Claudin 1 is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
Applications	IHC/IF
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SR0646
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC/IF (H, M, R) 1:500-1:1000
Formulation	PBS with 0.15% ProClin300, 100 Mu g/mL BSA and 50% glycerol.
Isotype	IgG
Storage Instruction	Store at -20C for up to one year, and avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9076
Gene Symbol	CLDN1
Uniprot ID	CLD1_HUMAN
Immunogen	KLH conjugated Synthetic peptide corresponding to Human claudin 1
Immunogen Region	
Specificity	
Immunogen Sequence	