

Anti-Recombinant-PD1 antibody [SM0656] (STJA0010656)

STJA0010656

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Recombinant-PD1 is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
Applications	IHC/IF
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SM0656
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC/IF (H) 1:500-1:1000
Formulation	PBS with 0.15% ProClin300, 100 µg/mL BSA and 50% glycerol.
Isotype	IgG1
Storage Instruction	Store at -20°C for up to one year, and avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Gene ID	18566
Gene Symbol	Pdcd1
Uniprot ID	PDCD1_MOUSE
Immunogen	KLH conjugated Synthetic peptide corresponding to Mouse PD1
Immunogen Region	
Specificity	
Immunogen Sequence	