

Anti-Recombinant-Bcl-2 antibody [SR0674] (STJA0010674)

STJA0010674

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit monoclonal antibody anti-Recombinant-Bcl-2 is suitable for use in Western Blot, Immunohistochemistry and Immunofluorescence research applications.
Applications	WB/IHC/IF
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SR0674
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB (H, M) 1:500-1:1000 IHC/IF (H, M) 1:500-1:1000
Formulation	PBS with 0.15% ProClin300, 100 Mu g/mL BSA and 50% glycerol.
Isotype	IgG
Storage	Store at -20C for up to one year, and avoid repeated freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	596
Gene Symbol	BCL2
Uniprot ID	BCL2_HUMAN
Immunogen	Recombinant protein corresponding to Human Bcl-2
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	