

Anti-Recombinant-RFP antibody [SM0679] (STJA0010679)

STJA0010679

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Recombinant-RFP is suitable for use in Immunocytochemistry and Immunofluorescence research applications.
Applications	ICC/IF
Host/Source	Mouse
Reactivity	RFP

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SM0679
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	ICC/IF (Species independent) 1:500-1:1000
Formulation	PBS with 0.15% ProClin300, 100 µg/mL BSA and 50% glycerol.
Isotype	IgG1
Storage Instruction	Store at -20°C for up to one year, and avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Recombinant protein corresponding to RFP
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	