

Anti-Recombinant-Aquaporin 4 antibody [SM0682] (STJA0010682)

STJA0010682

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Recombinant-Aquaporin 4 is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
Applications	IHC/IF
Host/Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SM0682
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC/IF (H, M, R) 1:250-1:1000
Formulation	PBS with 0.15% ProClin300, 100 µg/mL BSA and 50% glycerol.
Isotype	IgG1
Storage	Store at -20°C for up to one year, and avoid repeated freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	11829
Gene Symbol	Aqp4
Uniprot ID	AQP4_MOUSE
Immunogen	Recombinant protein corresponding to Mouse Aquaporin 4
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	