

Anti-Recombinant-EIF2S1 antibody [SM0684] (STJA0010684)

STJA0010684

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Recombinant-EIF2S1 is suitable for use in Western Blot, Immunohistochemistry and Immunofluorescence research applications.
Applications	WB/IHC/IF
Host/Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SM0684
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB (H, M, R) 1:500-1:1000 IHC/IF (M, R) 1:2000-1:8000
Formulation	PBS with 0.15% ProClin300, 100 µg/mL BSA and 50% glycerol.
Isotype	IgG1
Storage	Store at -20°C for up to one year, and avoid repeated freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	13665
Gene Symbol	Eif2s1
Uniprot ID	IF2A_MOUSE
Immunogen	Recombinant protein corresponding to Mouse EIF2S1
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	