

Anti-Recombinant-Bax antibody [SR0701] (STJA0010701)

STJA0010701

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit monoclonal antibody anti-Recombinant-Bax is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SR0701
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB (H, M, R) 1:1000-1:2000
Formulation	PBS with 0.15% ProClin300, 100 Mu g/mL BSA and 50% glycerol.
Isotype	IgG
Storage Instruction	Store at -20C for up to one year, and avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Gene ID	12028
Gene Symbol	Bax
Uniprot ID	BAX_MOUSE
Immunogen	KLH conjugated Synthetic peptide corresponding to Mouse Bax
Immunogen Region	
Specificity	
Immunogen Sequence	