

Anti-Recombinant-160 kDa Neurofilament Medium antibody [SM0702] (STJA0010702)

STJA0010702

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Recombinant-160 kDa Neurofilament Medium is suitable for use in Western Blot, Immunohistochemistry and Immunofluorescence research applications.
Applications	WB/IHC/IF
Host/Source	Mouse
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SM0702
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB (M, R) 1:500-1:2000 IHC/IF (H, M, R, Mk) 1:500-1:1000
Formulation	PBS with 0.15% ProClin300, 100 Mu g/mL BSA and 50% glycerol.
Isotype	IgG1
Storage	Store at -20C for up to one year, and avoid repeated freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	18040
Gene Symbol	Nefm
Uniprot ID	NFM_MOUSE
Immunogen	Recombinant protein corresponding to Mouse 160 kDa Neurofilament Medium
Immunogen Region	
Specificity	
Immunogen Sequence	