

Anti-Recombinant-Cytokeratin 6a antibody [SM0703] (STJA0010703)

STJA0010703

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Recombinant-Cytokeratin 6a is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
Applications	IHC/IF
Host/Source	Mouse
Reactivity	Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SM0703
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC/IF (M, R) 1:500-1:1000
Formulation	PBS with 0.15% ProClin300, 100 Mu g/mL BSA and 50% glycerol.
Isotype	IgG1
Storage Instruction	Store at -20C for up to one year, and avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Gene ID	16687
Gene Symbol	Krt6a
Uniprot ID	K2C6A_MOUSE
Immunogen	Recombinant protein corresponding to Mouse Cytokeratin 6a
Immunogen Region	
Specificity	
Immunogen Sequence	