

**Anti-Recombinant-Caldesmon/CDM antibody [SM0714] (STJA0010714)**

STJA0010714

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-Recombinant-Caldesmon/CDM is suitable for use in Western Blot, Immunohistochemistry and Immunofluorescence research applications.
<b>Applications</b>	WB/IHC/IF
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	SM0714
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity purification
<b>Dilution Range</b>	WB (H) 1:500-1:1000 IHC/IF (H) 1:1000-1:2000
<b>Formulation</b>	PBS with 0.15% ProClin300, 100 µg/mL BSA and 50% glycerol.
<b>Isotype</b>	IgG1
<b>Storage</b>	Store at -20°C for up to one year, and avoid repeated freeze-thaw cycles.
<b>Instruction</b>	

**TARGET INFORMATION**

<b>Gene ID</b>	<a href="#">800</a>
<b>Gene Symbol</b>	<a href="#">CALD1</a>
<b>Uniprot ID</b>	<a href="#">CALD1_HUMAN</a>
<b>Immunogen</b>	Recombinant protein corresponding to Human Caldesmon
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	