

Anti-Recombinant-MYC-Tag antibody [SM0719] {HRP} (STJA0010719)

STJA0010719

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Recombinant-MYC-Tag is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Mouse
Reactivity	Myc-Tag

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SM0719
Concentration	
Conjugation	HRP
Purification	Affinity purification
Dilution Range	WB (Species independent) 1:500-1:1000
Formulation	PBS with 0.15% ProClin300, 100 Mu g/mL BSA and 50% glycerol.
Isotype	
Storage Instruction	Store at -20C for up to one year, and avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	KLH conjugated Synthetic peptide corresponding to EQKLISEEDL
Immunogen Region	
Specificity	
Immunogen Sequence	