

Anti-Recombinant-His-Tag antibody [SR0722] {HRP} (STJA0010722)

STJA0010722

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit monoclonal antibody anti-Recombinant-His-Tag is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	His-Tag

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SR0722
Concentration	
Conjugation	HRP
Purification	Affinity purification
Dilution Range	WB (Species independent) 1:2000-1:5000
Formulation	PBS with 0.15% ProClin300, 100 Mu g/mL BSA and 50% glycerol.
Isotype	
Storage Instruction	Store at -20C for up to one year, and avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	KLH conjugated Synthetic peptide corresponding to HHHHHH (6* His)
Immunogen Region	
Specificity	
Immunogen Sequence	