

Anti-Recombinant-beta Tubulin antibody [SM0724] {HRP} (STJA0010724)

STJA0010724

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Recombinant-beta Tubulin is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SM0724
Concentration	
Conjugation	HRP
Purification	Affinity purification
Dilution Range	WB (H, M, R) 1:1000-1:2000
Formulation	PBS with 0.15% ProClin300, 100 Mu g/mL BSA and 50% glycerol.
Isotype	
Storage Instruction	Store at -20C for up to one year, and avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Gene ID	203068
Gene Symbol	TUBB
Uniprot ID	TBB5_HUMAN
Immunogen	Recombinant protein corresponding to Human beta Tubulin
Immunogen Region	
Specificity	
Immunogen Sequence	