

Anti-Recombinant-alpha Tubulin antibody [SM0725] {HRP} (STJA0010725)

STJA0010725

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Recombinant-alpha Tubulin is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SM0725
Concentration	
Conjugation	HRP
Purification	Affinity purification
Dilution Range	WB (H, M, R) 1:3000-1:6000
Formulation	PBS with 0.15% ProClin300, 100 Mu g/mL BSA and 50% glycerol.
Isotype	
Storage Instruction	Store at -20C for up to one year, and avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Gene ID	22142
Gene Symbol	Tuba1a
Uniprot ID	TBA1A_MOUSE
Immunogen	Recombinant protein corresponding to Mouse alpha Tubulin
Immunogen Region	
Specificity	
Immunogen Sequence	