

Anti-NR1D1 antibody (591-605 aa) [R05-8B1] (STJA0011699)

STJA0011699

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit monoclonal antibody anti-NR1D1 (591-605 aa) is suitable for use in Western Blot, Immunocytochemistry and Immunofluorescence research applications.
Applications	WB/ICC/IF
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Clone ID	R05-8B1
Concentration	0.3 mg/ml
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000 IF 1:50-1:200
Formulation	50mM Tris-Glycine (pH7.4) , 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Isotype	IgG
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	9572
Gene Symbol	NR1D1
Uniprot ID	NR1D1_HUMAN
Immunogen	A synthetic peptide of human NR1D1
Immunogen Region	591-605 aa
Specificity	
Immunogen Sequence	



Immunocytochemistry analysis of NR1D1 (green) in HeLa using NR1D1 antibody, and DAPI (blue)



Western blot analysis of NR1D1 in K562, C6, 3T3, HeLa lysates using NR1D1 antibody.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081