

Anti-BAG4 antibody (355-369 aa) [R06-6J7] (STJA0011961)

STJA0011961

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit monoclonal antibody anti-SODD (355-369 aa) is suitable for use in Western Blot, Immunohistochemistry, Immunocytochemistry and Immunofluorescence research applications.
Applications	WB/IHC-F/IHC-P/ICC/IF
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	R06-6J7
Concentration	0.3 mg/ml
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000 IHC 1:50-1:100 IF 1:50-1:200
Formulation	50mM Tris-Glycine (pH7.4) , 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Isotype	IgG
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	9530
Gene Symbol	BAG4
Uniprot ID	BAG4_HUMAN
Immunogen	A synthetic peptide of human Silencer of Death Domain
Immunogen Region	355-369 aa
Specificity	
Immunogen Sequence	



Western blot analysis of Silencer of Death Domain in MOLT4 lysates using SODD antibody.