

Anti-TIMM44 antibody (30-229 aa) [R04-7H8] (STJA0012056)

STJA0012056

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit monoclonal antibody anti-TIMM44 (30-229 aa) is suitable for use in Western Blot, Immunocytochemistry, Immunofluorescence and Immunoprecipitation research applications.
Applications	WB/ICC/IF/IP
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Clone ID	R04-7H8
Concentration	0.3 mg/ml
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000 IF 1:50-1:200 IP 1:20
Formulation	50mM Tris-Glycine (pH7.4) , 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Isotype	IgG
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION

Gene ID	10469
Gene Symbol	TIMM44
Uniprot ID	TIM44_HUMAN
Immunogen	Recombinant protein of human TIMM44
Immunogen Region	30-229 aa
Specificity	
Immunogen Sequence	



Immunocytochemistry analysis of TIMM44 (green) in HeLa using TIMM44 antibody, and DAPI (blue)



Western blot analysis of TIMM44 in HeLa lysates using TIMM44 antibody.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081