

Anti-XDH antibody (1200-1300 aa) [R07-1D5] (STJA0012141)

STJA0012141

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit monoclonal antibody anti-Xanthine Oxidase (1200-1300 aa) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	R07-1D5
Concentration	0.3 mg/ml
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000
Formulation	50mM Tris-Glycine (pH7.4) , 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Isotype	IgG
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION

Gene ID	7498
Gene Symbol	XDH
Uniprot ID	XDH_HUMAN
Immunogen	A synthetic peptide of human Xanthine Oxidase
Immunogen Region	1200-1300 aa
Specificity	
Immunogen Sequence	



Western blot analysis of Xanthine Oxidase in K562 lysates using Xanthine Oxidase antibody.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081