

Anti-PPP3CA antibody (1-14 aa) [R01-6F5] (STJA0012295)

STJA0012295

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit monoclonal antibody anti-Calcineurin A (1-14 aa) is suitable for use in Western Blot and Immunoprecipitation research applications.
Applications	WB/IP
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Clone ID	R01-6F5
Concentration	0.3 mg/ml
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000 IP 1:20
Formulation	50mM Tris-Glycine (pH7.4) , 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Isotype	IgG
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION

Gene ID	5530
Gene Symbol	PPP3CA
Uniprot ID	PP2BA_HUMAN
Immunogen	A synthetic peptide of human Calcineurin A
Immunogen Region	1-14 aa
Specificity	
Immunogen Sequence	



Western blot analysis of Calcineurin A in rat brain, mouse lung lysates using Calcineurin A antibody.



Western blot analysis of Calcineurin A in HeLa, A549, U2OS lysates using Calcineurin A antibody.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081