

Anti-LGALS1 antibody (97-110 aa) [R04-0C8] (STJA0012383)

STJA0012383

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit monoclonal antibody anti-Galectin 1 (97-110 aa) is suitable for use in Western Blot and Immunoprecipitation research applications.
Applications	WB/IP
Host/Source	Rabbit
Reactivity	Mouse/Rat/Hamster

PRODUCT PROPERTIES


Clonality	Monoclonal
Clone ID	R04-0C8
Concentration	0.3 mg/ml
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000 IP 1:20
Formulation	50mM Tris-Glycine (pH7.4) , 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Isotype	IgG
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION

Gene ID	3956
Gene Symbol	LGALS1
Uniprot ID	LEG1_HUMAN
Immunogen	A synthetic peptide of human Galectin 1
Immunogen Region	97-110 aa
Specificity	
Immunogen Sequence	



Western blot analysis of Galectin 1 in 3T3 lysates using Galectin 1 antibody.



Western blot analysis of Galectin 1 in CHO-K1, C6 lysates using Galectin 1 antibody.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081