

Anti-GRK6 antibody (97-110 aa) [R07-3E4] (STJA0012396)

STJA0012396

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit monoclonal antibody anti-GRK6 (97-110 aa) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Clone ID	R07-3E4
Concentration	0.3 mg/ml
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000
Formulation	50mM Tris-Glycine (pH7.4) , 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Isotype	IgG
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION

Gene ID	2870
Gene Symbol	GRK6
Uniprot ID	GRK6_HUMAN
Immunogen	A synthetic peptide of human GRK6
Immunogen Region	97-110 aa
Specificity	
Immunogen Sequence	



Western blot analysis of GRK6 in mouse spleen lysates using GRK6 antibody.



Western blot analysis of GRK6 in A549, HL-60, U2OS lysates using GRK6 antibody.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081