

Anti-THBS1 antibody (547-1170 aa) [R06-5F3] (STJA0012618)

STJA0012618

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit monoclonal antibody anti-Thrombospondin 1 (547-1170 aa) is suitable for use in Western Blot, Immunohistochemistry and Immunoprecipitation research applications.
Applications	WB/IHC-P/IP
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES


Clonality	Monoclonal
Clone ID	R06-5F3
Concentration	0.3 mg/ml
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000 IHC 1:50-1:100 IP 1:20
Formulation	50mM Tris-Glycine (pH7.4) , 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Isotype	IgG
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	7057
Gene Symbol	THBS1
Uniprot ID	TSP1_HUMAN
Immunogen	Recombinant protein of human Thrombospondin 1
Immunogen Region	547-1170 aa
Specificity	
Immunogen Sequence	



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using Thrombospondin 1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of Thrombospondin 1 in MCF-7, 293T, C2C12 lysates using Thrombospondin 1 antibody.