

## Anti-CD81 antibody (34-63+113-201 aa) [R02-3D6] (STJA0013040)

STJA0013040

### GENERAL INFORMATION


<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit monoclonal antibody anti-CD81 (34-63+113-201 aa) is suitable for use in Western Blot, Immunohistochemistry, Immunocytochemistry and Immunofluorescence research applications.
<b>Applications</b>	WB/IHC-F/IHC-P/ICC/IF
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human/Mouse

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	R02-3D6
<b>Concentration</b>	0.3 mg/ml
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity Purified
<b>Dilution Range</b>	WB 1:500-1:1000 IHC 1:50-1:100 IF 1:50-1:200
<b>Formulation</b>	50mM Tris-Glycine (pH7.4) , 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
<b>Isotype</b>	IgG
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	975
<b>Gene Symbol</b>	CD81
<b>Uniprot ID</b>	CD81_HUMAN
<b>Immunogen</b>	Recombinant protein of human CD81
<b>Immunogen Region</b>	34-63+113-201 aa
<b>Specificity</b>	
<b>Immunogen Sequence</b>	



Western blot analysis of CD81 in C2C12, 3T3, MCF-7 lysates using CD81 antibody.