

Anti-HSF1 antibody (485-529 aa) [R09-3M8] (STJA0014242)

STJA0014242

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit monoclonal antibody anti-Hsf1 (485-529 aa) is suitable for use in Western Blot, Immunohistochemistry, Immunocytochemistry, Immunofluorescence, Immunoprecipitation and Flow Cytometry research applications.
Applications	WB/IHC-P/ICC/IF/IP/FC
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	R09-3M8
Concentration	-
Conjugation	Unconjugated
Purification	Affinity Chromatography
Dilution	WB 1:500-1:1000
Range	IHC 1:50-1:100 IF 1:50-1:200 IP 1:50 FC 1:50-1:100
Formulation	Liquid in PBS pH7.4, 150mM NaCl, 0.02% Sodium Azide and 50% Glycerol.
Isotype	IgG
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	3297
Gene Symbol	HSF1
Uniprot ID	HSF1_HUMAN
Immunogen	A synthesized peptide derived from human HSF1
Immunogen Region	485-529 aa
Specificity	
Immunogen Sequence	



Western blot analysis of HSF1 in HeLa lysates using Hsf1 antibody.