

Anti-TRADD antibody (210-312 aa) [R03-7G5] (STJA0014295)

STJA0014295

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit monoclonal antibody anti-TRADD (210-312 aa) is suitable for use in Western Blot, Immunohistochemistry, Immunocytochemistry, Immunofluorescence and Immunoprecipitation research applications.
Applications	WB/IHC-P/ICC/IF/IP
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	R03-7G5
Concentration	-
Conjugation	Unconjugated
Purification	Affinity Chromatography
Dilution	WB 1:500-1:1000
Range	IHC 1:50-1:100 IF 1:50-1:200 IP 1:10
Formulation	Liquid in PBS pH7.4, 150mM NaCl, 0.02% Sodium Azide and 50% Glycerol.
Isotype	IgG
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION

Gene ID	8717
Gene Symbol	TRADD
Uniprot ID	TRADD_HUMAN
Immunogen	A synthesized peptide derived from human TRADD
Immunogen Region	210-312 aa
Specificity	
Immunogen Sequence	