

Anti-MUC4 antibody (1564-1673 aa) [R02-7G1] (STJA0014931)

STJA0014931

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit monoclonal antibody anti-Mucin 4 (1564-1673 aa) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	R02-7G1
Concentration	-
Conjugation	Unconjugated
Purification	Affinity Chromatography
Dilution Range	WB 1:500-1:1000
Formulation	Liquid in PBS pH7.4, 150mM NaCl, 0.02% Sodium Azide and 50% Glycerol.
Isotype	IgG
Storage Instruction	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.

TARGET INFORMATION

Gene ID	4585
Gene Symbol	MUC4
Uniprot ID	MUC4_HUMAN
Immunogen	A synthesized peptide derived from human MUC4
Immunogen Region	1564-1673 aa
Specificity	
Immunogen Sequence	



Western blot analysis of MUC4 in MCF-7 lysates using Mucin 4 antibody.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081