

Anti-LGI1 antibody (150-250 aa) [R03-8B6] (STJA0015620)

STJA0015620

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit monoclonal antibody anti-LGI1 (150-250 aa) is suitable for use in Western Blot, Immunohistochemistry and Immunoprecipitation research applications.
Applications	WB/IHC-P/IP
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	R03-8B6
Concentration	-
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:500-1:1000 IHC 1:50-1:100 IP 1:50
Formulation	Liquid in PBS pH7.4, 150mM NaCl, 0.02% Sodium Azide and 50% Glycerol.
Isotype	IgG
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	9211
Gene Symbol	LGI1
Uniprot ID	LGI1_HUMAN
Immunogen	A synthesized peptide derived from human Lgi1
Immunogen Region	150-250 aa
Specificity	
Immunogen Sequence	



Western blot analysis of Lgi1 in HeLa lysates using LGI1 antibody.