

Anti-PEX19 antibody (30-89 aa) [R02-4Y2] (STJA0015623)

STJA0015623

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit monoclonal antibody anti-PEX19 (30-89 aa) is suitable for use in Western Blot, Immunocytochemistry, Immunofluorescence, Immunoprecipitation and Flow Cytometry research applications.
Applications	WB/ICC/IF/IP/FC
Host/Source	Rabbit
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	R02-4Y2
Concentration	-
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution	WB 1:500-1:1000
Range	IF 1:50-1:200 IP 1:50 FC 1:50-1:100
Formulation	Liquid in PBS pH7.4, 150mM NaCl, 0.02% Sodium Azide and 50% Glycerol.
Isotype	IgG
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	5824
Gene Symbol	PEX19
Uniprot ID	PEX19_HUMAN
Immunogen	A synthesized peptide derived from human PEX19
Immunogen Region	30-89 aa
Specificity	
Immunogen Sequence	



Western blot analysis of PEX19 in Jurkat lysates using PEX19 antibody.