

Anti-CD329 antibody [K8] {PE/Cyanine 7} (STJA0016234)

STJA0016234

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Sialic Acid-Binding Ig-Like Lectin 9 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	K8
Concentration	5 Mu L/Tests
Conjugation	PE/Cyanine 7
Purification	Affinity purified
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

TARGET INFORMATION

Gene ID	27180
Gene Symbol	SIGLEC9
Uniprot ID	SIGL9_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	