

Anti-CD371 antibody [50C1] {FITC} (STJA0016238)

STJA0016238

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-C-Type Lectin Domain Family 12 Member A is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	50C1
Concentration	5 Mu L/Tests
Conjugation	FITC
Purification	Affinity purified
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

TARGET INFORMATION

Gene ID	160364
Gene Symbol	CLEC12A
Uniprot ID	CL12A_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	