

**Anti-CD41 antibody [10E5] {PE/Cyanine 7} (STJA0016251)**

STJA0016251

**GENERAL INFORMATION**

|                          |                                                                                                                 |
|--------------------------|-----------------------------------------------------------------------------------------------------------------|
| <b>Product Type</b>      | Primary antibodies                                                                                              |
| <b>Short Description</b> | Mouse recombinant antibody anti-Integrin Alpha-Iib is suitable for use in Flow Cytometry research applications. |
| <b>Applications</b>      | FC                                                                                                              |
| <b>Host/Source</b>       | Mouse                                                                                                           |
| <b>Reactivity</b>        | Human                                                                                                           |

**PRODUCT PROPERTIES**

|                            |                                                                                                                         |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------|
| <b>Clonality</b>           | Recombinant                                                                                                             |
| <b>Clone ID</b>            | 10E5                                                                                                                    |
| <b>Concentration</b>       | 5 µg L/Tests                                                                                                            |
| <b>Conjugation</b>         | PE/Cyanine 7                                                                                                            |
| <b>Purification</b>        | Affinity purified                                                                                                       |
| <b>Dilution Range</b>      |                                                                                                                         |
| <b>Formulation</b>         | Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.                             |
| <b>Isotype</b>             | IgG2ak                                                                                                                  |
| <b>Storage Instruction</b> | This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze. |

**TARGET INFORMATION**

|                           |                             |
|---------------------------|-----------------------------|
| <b>Gene ID</b>            | <a href="#">3674</a>        |
| <b>Gene Symbol</b>        | <a href="#">ITGA2B</a>      |
| <b>Uniprot ID</b>         | <a href="#">ITA2B_HUMAN</a> |
| <b>Immunogen</b>          |                             |
| <b>Immunogen Region</b>   |                             |
| <b>Specificity</b>        |                             |
| <b>Immunogen Sequence</b> |                             |