

Anti-CD41 antibody [10E5] {PerCP/Cyanine 5.5} (STJA0016252)

STJA0016252

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse recombinant antibody anti-Integrin Alpha-Iib is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Recombinant
Clone ID	10E5
Concentration	5 µg L/Tests
Conjugation	PerCP/Cyanine 5.5
Purification	Affinity purified
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

TARGET INFORMATION

Gene ID	3674
Gene Symbol	ITGA2B
Uniprot ID	ITA2B_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	