

Anti-CD44 antibody [HI44a] {PE/Cyanine 5} (STJA0016260)

STJA0016260

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Cd44 antigen is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	HI44a
Concentration	5 µg L/Tests
Conjugation	PE/Cyanine 5
Purification	Affinity purified
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2ak
Storage Instruction	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

TARGET INFORMATION

Gene ID	960
Gene Symbol	CD44
Uniprot ID	CD44_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	