

Anti-IgG2aL antibody [G013C12] {PE/EF594} (STJA0016270)
STJA0016270**GENERAL INFORMATION**

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-IgG2aL is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	G013C12
Concentration	5 Mu L/Test
Conjugation	PE/EF594
Purification	Affinity purified
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2aL
Storage Instruction	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen Sequence