

Anti-IgG2aL antibody [B39-4] {EF700} (STJA0016272)

STJA0016272

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-IgG2aL is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	B39-4
Concentration	
Conjugation	EF700
Purification	Affinity purified
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2aL
Storage Instruction	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	