

Anti-IgG antibody [SHG-1] {EV500} (STJA0016279)

STJA0016279

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Syrian Hamster monoclonal antibody anti-IgG is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Syrian Hamster
Reactivity	

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SHG-1
Concentration	
Conjugation	EV500
Purification	Affinity purified
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG
Storage Instruction	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen Sequence