

Anti-IgG2ak antibody [C1.18.4] {EV500} (STJA0016291)

STJA0016291

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Mouse monoclonal antibody anti-IgG2ak is suitable for use in Flow Cytometry research applications. |
| Applications | FC |
| Host/Source | Mouse |
| Reactivity | |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Monoclonal |
| Clone ID | C1.18.4 |
| Concentration | |
| Conjugation | EV500 |
| Purification | Affinity purified |
| Dilution Range | |
| Formulation | Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant. |
| Isotype | IgG2ak |
| Storage Instruction | This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze. |

TARGET INFORMATION

| |
|---------------------------|
| Gene ID |
| Gene Symbol |
| Uniprot ID |
| Immunogen |
| Immunogen Region |
| Specificity |
| Immunogen Sequence |