

Anti-MHC II antibody [M5/114] {EV500} (STJA0016312)

STJA0016312

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-MHC II is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	M5/114
Concentration	
Conjugation	EV500
Purification	Affinity purified
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

TARGET INFORMATION

Gene ID	14961
Gene Symbol	H2-Ab1
Uniprot ID	HB2A_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	