

Anti-CD11c antibody [N418] {EV500} (STJA0016315)

STJA0016315

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Armenian Hamster monoclonal antibody anti-Integrin Alpha-X is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Armenian Hamster
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	N418
Concentration	
Conjugation	EV500
Purification	Affinity purified
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG
Storage Instruction	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

TARGET INFORMATION

Gene ID	16411
Gene Symbol	Itgax
Uniprot ID	ITAX_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	