

Anti-CD3 antibody [17A2] {EV610} (STJA0016329)

STJA0016329

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rat monoclonal antibody anti-CD3 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Rat
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	17A2
Concentration	
Conjugation	EV610
Purification	Affinity purified
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

TARGET INFORMATION

Gene ID	12501 12502 12503 Cd3e Cd3g CD3E_MOUSE CD3G_MOUSE CD3Z_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	