

Anti-CD22 antibody [Cy34.1] {EF700} (STJA0016332)

STJA0016332

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-B-Cell Receptor Cd22 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	Cy34.1
Concentration	
Conjugation	EF700
Purification	Affinity purified
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

TARGET INFORMATION

Gene ID	12483
Gene Symbol	Cd22
Uniprot ID	CD22_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	