

Anti-CD32 antibody [IV-3] {PE/Cyanine 5.5} (STJA0016361)

STJA0016361

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Low Affinity Immunoglobulin Gamma Fc Region Receptor II-B is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	IV-3
Concentration	5 Mu L/Test
Conjugation	PE/Cyanine 5.5
Purification	Affinity purified
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG2bk
Storage Instruction	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

TARGET INFORMATION

Gene ID	2213
Gene Symbol	FCGR2B
Uniprot ID	FCG2B_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	