

Anti-CD15/SSEA-1 antibody [HI98] {EF700} (STJA0016365)

STJA0016365

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Alpha (1-3-Fucosyltransferase 4 is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	HI98
Concentration	
Conjugation	EF700
Purification	Affinity purified
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgMk
Storage Instruction	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

TARGET INFORMATION

Gene ID	2526
Gene Symbol	FUT4
Uniprot ID	FUT4_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen	
Sequence	