

Anti-CD64 antibody [10.1] {EF780} (STJA0016370)

STJA0016370

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-High Affinity Immunoglobulin Gamma Fc Receptor I is suitable for use in Flow Cytometry research applications.
Applications	FC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	10.1
Concentration	5 Mu L/Test
Conjugation	EF780
Purification	Affinity purified
Dilution Range	
Formulation	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
Isotype	IgG1k
Storage Instruction	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

TARGET INFORMATION

Gene ID	2209
Gene Symbol	FCGR1A
Uniprot ID	FCGR1_HUMAN
Immunogen	
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	