

**Anti-CD8a antibody [53-6.7] {EV500} (STJA0016388)**

STJA0016388

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rat monoclonal antibody anti-T-Cell Surface Glycoprotein Cd8 Alpha Chain is suitable for use in Flow Cytometry research applications.
<b>Applications</b>	FC
<b>Host/Source</b>	Rat
<b>Reactivity</b>	Mouse

**PRODUCT PROPERTIES**

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	53-6.7
<b>Concentration</b>	
<b>Conjugation</b>	EV500
<b>Purification</b>	Affinity purified
<b>Dilution Range</b>	
<b>Formulation</b>	Phosphate buffered solution, pH 7.2, containing 0.09% stabilizer and 1% protein protectant.
<b>Isotype</b>	IgG2ak
<b>Storage Instruction</b>	This product can be stored at 2-8°C for 12 months. Please protected from prolonged exposure to light and do not freeze.

**TARGET INFORMATION**

<b>Gene ID</b>	<a href="#">12525</a>
<b>Gene Symbol</b>	<a href="#">Cd8a</a>
<b>Uniprot ID</b>	<a href="#">CD8A_MOUSE</a>
<b>Immunogen</b>	
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen</b>	
<b>Sequence</b>	